



MEASUREMENT AND SENSOR
TECHNOLOGY



Department of
Electrical Engineering and
Information Technology



CHEMNITZ UNIVERSITY OF
TECHNOLOGY

Graphical user interface (GUI) on PC for communication with potentiostat embedded system

Rahul Chowdary Pinnaka (621023)

Project Lab Embedded Systems

Chemnitz, 13.07.2020

Outline

- Motivation
- Introduction
- Graphical user interface (GUI) for electro chemistry techniques
- Embedded platform
- Conclusion
- References

Motivation

- Reduce the manual effort to record the measurement readings.
- Increase the efficiency of the individual during the measurement process.
- Design a graphical user interface (GUI) to automate the measurement process
 - Record the measurement data
 - Plot real time graphs
 - Save the measurement data to excel file



Fig 1: Image showing the manual work during and after the experiment [7]

Introduction

- Design the GUI on personal computer (PC) for electro chemical techniques
 - Cyclic voltammetry
 - Square wave voltammetry and
 - Impedance spectroscopy
- Read out the measurement data from the potentiostat over Universal Asynchronous Receiver Transmitter (UART) communication protocol
- Plot the real time graph and store the readings in an excel file




Fig 2: Image showing the overview of the project

Graphical user interface (GUI)

- Based on Palmsens PSTrace software for electro chemical processes.
- Functions of PSTrace software:
 - Automate the electro chemical measurement process
 - Plotting the real time graphs of the measurement data
 - Analysis of the measurement data
 - Saving the data to excel file
- GUI consists of 4 pages namely:
 - Home page
 - Cyclic voltammetry
 - Square wave voltammetry
 - Impedance spectroscopy

Home page

- Contains brief introduction about the project and various features included in the GUI.




Fig 3: Home page of the GUI

Cyclic voltammetry

- Cyclic voltammetry is a type of electrochemical measurement technique.

TECHNISCHE UNIVERSITÄT CHEMNITZ
MST
MEASUREMENT AND SENSOR TECHNOLOGY

Electro-chemistry Techniques

Cyclic voltammetry

Square wave voltammetry

Impedance spectroscopy

Cyclic Voltammetry

E start V

E stop V

Current

E step V

Scan rate V/s

No of scans

START

Version 1.0

Copyright © 2020 TU Chemnitz MST Department - All rights reserved

Fig 4: Home page of cyclic voltammetry

Cyclic voltammetry


- Accepted values for the parameters are specified in the Table I

No	Parameter	Value
1	E_{start}	Between -1.5 V and +1.5 V
2	E_{stop}	Between -1.5 V and +1.5 V and should be greater than E_{start}
3	Current	Between 100 pA to 100 mA
4	E_{step}	Between 0.01 V and 0.1 V
5	Scan rate	Between 0.01 V and 0.1 V
6	No of scans	Between 1 and 5

Table I: Range specification for parameters in cyclic voltammetry

Cyclic voltammetry

- Data validations are carried out to make sure that no fields are left empty and data entered in the corresponding parameter fields is within the specified range.



The screenshot displays the 'Cyclic Voltammetry' software interface. On the left is a sidebar with the logo of Technische Universität Chemnitz and MST (Measurement and Sensor Technology). Below the logo, it lists 'Electro-chemistry Techniques' with three options: Cyclic voltammetry (selected), Square wave voltammetry, and Impedance spectroscopy. The main area is titled 'Cyclic Voltammetry' and contains several input fields with validation messages:

Parameter	Value	Unit	Validation Message
E start	2	V	* Enter E start between -1.5 V and +1.5 V
E stop	-3.5	V	* Enter E stop between -1.5 V and +1.5 V
Current			* Please select a value
E step		V	* Please select a value
Scan rate		V/s	* Please select a value
No of scans			* Please select a value

Below the input fields is a blue 'START' button. At the bottom of the interface, there is a footer bar with 'Version 1.0' on the left and 'Copyright © 2020 TU Chemnitz MST Department - All rights reserved' on the right.

Fig 5: Range specific validations

Cyclic voltammetry

- Parameters can be sent to the microcontroller when the data entered in the fields is within in range specified and the experiment can be started.

The screenshot shows a software window titled "CIFE". On the left is a sidebar with the logo of Technische Universität Chemnitz and MST (Measurement and Sensor Technology). Below the logo is a blue bar with the text "Electro-chemistry Techniques". Underneath are three options: "Cyclic voltammetry" (selected), "Square wave voltammetry", and "Impedance spectroscopy". At the bottom of the sidebar is a blue bar with "Version 1.0". The main area is titled "Cyclic Voltammetry" and contains several input fields: "E start" with the value "-1" and unit "V", "E stop" with the value "1" and unit "V", "Current" with a dropdown menu showing "10 μ A", "E step" with a dropdown menu showing "0.05" and unit "V", "Scan rate" with a dropdown menu showing "0.01" and unit "V/s", and "No of scans" with a dropdown menu showing "1". Below these fields is a blue "START" button. At the bottom of the window is a blue bar with the text "Copyright © 2020 TU Chemnitz MST Department - All rights reserved".

Fig 6: Correct data entered in the parameter fields

Measurement data generation

- For cyclic voltammetry and square wave voltammetry techniques, the measurement current is generated using Ohm's law.
- Constant resistance of $10\text{k}\Omega$ is used in the measurement data generation.
- Initial value of measured voltage is the value of the parameter E_{start} and it increases until E_{stop} in steps of E_{step} value.
- Measurement current is calculated for each value of measurement voltage using ohm's law.
- For impedance spectroscopy technique, the measurement data of frequency and impedance is randomly generated.




Fig 7: Sample graph in cyclic voltammetry method

Cyclic voltammetry

- GUI receives output potential measured between the working electrode and the reference electrode, and the output current measured between the working electrode and the counter electrode as measurement data and plots the real time graph.




Fig 8: Real time graph of the received data

Cyclic voltammetry

- Real time graph can be saved as an image and measurement data can be exported into an excel file.




Fig 9: Real time graph saved as an image file

	A	B	C	D
1	Voltage (V)	Current (µA)		
2	-1.00004	-46.3724		
3	-0.999262	-45.6095		
4	-0.99848	-44.4174		
5	-0.997699	-43.8452		
6	-0.996996	-42.9392		
7	-0.996215	-42.9512		
8	-0.995434	-42.5936		
9	-0.994652	-41.6304		
10	-0.993949	-41.1393		

Fig 10: Real time graph data saved into an excel file

Square wave voltammetry

- Square wave voltammetry is a form of pulse voltammetry method available in electro chemistry techniques.
- Input potential is a combination of a square wave and staircase potential applied to the working electrode.
- The current at the working electrode is plotted versus the applied voltage to give the square wave voltammogram trace.

TECHNISCHE UNIVERSITÄT
CHEMNITZ
MST
MEASUREMENT AND SENSOR
TECHNOLOGY

**Electro-chemistry
Techniques**

Cyclic voltammetry

Square wave voltammetry

Impedance spectroscopy

Version 1.0

Square Wave Voltammetry

E start V

E stop V

Current ▾

E step ▾ V

Amplitude V

Frequency ▾ Hz

START

Copyright © 2020 TU Chemnitz MST Department - All rights reserved

Fig 11: Home page of square wave voltammetry

Square wave voltammetry

- Accepted values for the parameters are specified in the Table II.

No	Parameter	Value
1	E_{start}	Between -1.5 V and +1.5 V
2	E_{stop}	Between -1.5 V and +1.5 V and should be greater than E_{start}
3	Current	Between 100 pA to 100 mA
4	E_{step}	Between 0.01 V to 1 V
5	Amplitude	Between 0.1 V to 1 V
6	Frequency	Between 1 Hz and 150 Hz

Table II: Range specification for parameters

Impedance spectroscopy

- Electrochemical Impedance Spectroscopy is one of the most complex techniques available in electrochemical research.
- This technique is widely used to measure the impedance of a system in dependence of the AC potentials frequency.

TECHNISCHE UNIVERSITÄT CHEMNITZ
MST
MEASUREMENT AND SENSOR TECHNOLOGY

Electro-chemistry Techniques

- Cyclic voltammetry
- Square wave voltammetry
- Impedance spectroscopy**

Version 1.0

Impedance Spectroscopy

Scan type: Fixed

E dc: 0 V

E ac: mV

Current: V

Frequency type: V

Frequency: Hz

Min. frequency: V

Max. frequency: V

START

Copyright © 2020 TU Chemnitz MST Department - All rights reserved

Fig 12: Home page of Impedance spectroscopy

Impedance spectroscopy

- Accepted values for the parameters are specified in the Table III.

No	Parameter	Value
1	Scan type	Fixed
2	E_{dc}	0 V
3	E_{ac}	Between 10 mV and 50 mV
4	Current	Between 100 pA to 100 mA
5	Frequency type	Fixed (or) Scan
6	Frequency	Between 1 Hz and 150 Hz
7	Min frequency	Between 0.1 Hz and 100 kHz
8	Max frequency	100 kHz

Table III: Range specification for parameters

Embedded platform

- Hardware used in this project is STM32 Nucleo-L4R5ZI board.
- Integrated development environment (IDE) used for configuration of peripherals, generation of code, compilation and debugging is STM32 CubeIDE.
- The peripherals Timer and Universal Asynchronous Receiver Transmitter (UART) are used in this project.
- LPUART1 (Lower Power UART) is used in this project with Baud Rate set to 9600bits/sec, words length set to 8 bits




Fig 13: STM32 L4R5ZI board layout [6]

Conclusion

- GUI was designed successfully to select different electro chemistry methods and adjusted their input parameters and sent them to the microcontroller.
- Microcontroller board was configured and programmed to read the parameters from the GUI and sent the measurement data to the GUI.
- Input parameters in the GUI were non editable while receiving the measurement data.
- GUI received the measurement data and correctly interpreted the data and simultaneously plotted the real time graph.
- Features to save the graph as an image file and to export the graph data into an excel file were implemented.
- Extensive debugging and unit testing of the software was carried out to ensure the smooth operation of the application.
- Easy and versatile GUI was developed with advanced features suitable for all levels of user experience.

References

- [1] STMicroelectronics. Discovery Kit with STM32 Nucleo-L4R5ZI MCU User Manual. April 2018. UM1570.
- [2] STM32CubeIDE for STM32 configuration and initialization C code generation. March 2018. UM1718.
- [3] Getting started with STM32 MCU Discovery Kits software development tools. May 2016. UM2052.
- [4] Description of STM32 HAL and low-layer drivers User Manual. July 2017. UM1786.
- [5] Norris, Donald J. Programming with STM32 Getting started with Nucleo Board and C/C++. s.l.: McGraw-Hill Education, 2018.
- [6] Docs.zephyrproject.org. 2020. ST Nucleo L4R5ZI — Zephyr Project Documentation.
[online] Available at:
<https://docs.zephyrproject.org/latest/boards/arm/nucleo_l4r5zi/doc/index.html>
[Accessed 28 June 2020].
- [7] Bemidjistate.edu. 2020. Facilities & Equipment | Chemistry | Bemidji State University.
[online] Available at:
<<https://www.bemidjistate.edu/academics/departments/chemistry/about/facilities-equipment/>> [Accessed 4 July 2020].

Thank you



Group 11

Chair of Measurement and Sensor Technology



CHEMNITZ UNIVERSITY
OF TECHNOLOGY

20/20